

# Terahertz components

## Fiber-coupled THz emitter and detector



### FEATURES

- Based on unique GaBiAs photoconductive material
- Optimized for wavelengths around 1060 nm
- Pump pulse fiber delivery
- Technical passport and test report included

Fiber-coupled THz emitters and detectors offered by EKSPLA are designed for broadband operation and can be used in standard time-domain (THz-TDS) setup (Fig. 1). Unique GaBiAs

THz emitters and detectors are mounted into compact housing compatible with standard 1" optical holders. Performance of each device is checked and technical passport, including testing report, is provided for customer.




Fig. 1. THz time-domain spectroscopy setup

material used in photoconductive antennas (PCA) features excellent sensitivity for wavelengths up to approx. 1060 nm and electron lifetime shorter than 1 ps. As a result ultra short half-cycle THz pulses with broad spectra up to 5 THz can be generated and detected. Pump beam fiber delivery eliminates time-consuming adjustments and ensures maximum flexibility of experiment. As an example, this feature allows fast and convenient switching between different geometries: transmission, reflection, etc.

### INTEGRATED SILICON LENS

THz emitters and detectors are supplied with integrated hyper-hemispherical lenses, made from high-resistivity silicon, attached to PCA to increase the radiation efficiency of THz waves into free space. EKSPLA offers two standard types of these lenses: for collimated or diverging THz beam output. Advantage of collimated THz beam output is simple setup, because no additional optical components between THz emitter and detector are required for experiment. However, this design features bigger aberrations of THz beam, which affects focusing. In second case design of lens assumes positioning of the PCA in aplanatic point, which significantly reduces aberrations. As a result nearly diffraction limited spot of THz beam can be achieved.

### APPLICATIONS

- Time-resolved broadband THz spectroscopy
- Optical pump-THz probe spectroscopy
- Suitable for all-in-fiber system solution

## PHOTOCONDUCTIVE ANTENNA (PCA)

Photoconductive antennas are particularly designed for THz emitter or THz detector. The substrate of GaAs contains mesa-etched epitaxial active layer of GaBiAs in order to achieve high dark resistance. High photosensitivity of the material allows use of low average power optical pulses for excitation. On its surface a coplanar Hertzian type dipole antenna structure is formed using AuGeNi metallization (Fig. 2). The gap between metallic contacts is similar to laser spot diameter in detector case and larger – in emitter case. Photoconductive chip is mounted on PCB and placed inside metallic housing of device. SMA sockets on back side of the housing are used to connect DC or AC bias to THz emitter or lock-in amplifier input to THz detector.




Fig. 2. Microstip antenna drawings: (a) emitter, (b) detector (all dimensions are in micrometers)

## SPECIFICATIONS

| MODEL                                         | EMITTER EMT-10-FC                                 | DETECTOR DET-10-FC |
|-----------------------------------------------|---------------------------------------------------|--------------------|
| <b>PHOTOCONDUCTIVE ANTENNA</b>                |                                                   |                    |
| Photoconductive material                      | LT-GaBiAs                                         |                    |
| Dimensions of the wafer                       | 5 × 1.5 mm                                        |                    |
| Thickness                                     | 600 μm                                            |                    |
| Antenna type                                  | strip line                                        | dipole             |
| Photosensitivity                              | up to 1100 nm                                     |                    |
| Bias voltage                                  | 50 V max, 40 V typical                            | –                  |
| Detected THz bandwidth <sup>1)</sup>          | –                                                 | >4 THz             |
| <b>INTEGRATED FOCUSING LENS</b>               |                                                   |                    |
| Material                                      | HRFZ-silicon                                      |                    |
| Geometrical form                              | hyper-hemi-sphere                                 |                    |
| THz beam output                               | collimated or diverging                           | –                  |
| <b>FIBER DELIVERY</b>                         |                                                   |                    |
| Fiber length                                  | 3±0.2 m                                           |                    |
| Fiber connector                               | FC/APC                                            |                    |
| Wavelength <sup>2)</sup>                      | 1064 nm                                           |                    |
| Maximum optical pulse energy (on fiber input) | 1 nJ (30 mW at 30 MHz)                            |                    |
| Recommended pump source                       | EKSPLA LightWire FF50 (with double output option) |                    |

<sup>1)</sup> Pumped by 130 fs, 20 mW, 30 MHz pulses

<sup>2)</sup> Other wavelengths are available on request

## ORDERING INFORMATION

| DESCRIPTION                                       | MODEL     | NOTES                                                                                      |
|---------------------------------------------------|-----------|--------------------------------------------------------------------------------------------|
| Fiber-coupled THz emitter for 1060 nm wavelength  | EMT-10-FC | Includes Si lens, optical fiber with FC/APC connector and coaxial cable with BNC connector |
| Fiber-coupled THz detector for 1060 nm wavelength | DET-10-FC | Includes Si lens, optical fiber with FC/APC connector and coaxial cable with BNC connector |

## FIBER DELIVERY

Each fiber-coupled THz emitter and detector is delivered with specially design optical fiber. Standard configuration requires pre-chirped pulse with 0.5-2 ps duration on fiber input. Typically such pulse parameters are available directly from fiber lasers before compressor stage. In this case pulse is compressed while propagating through optical fiber with negative dispersion. Another solution can be applied for femtosecond pulses – near zero dispersion optical fiber. In such fiber femtosecond pulse keeps its form and pulse duration doesn't differ much on input and output of fiber. In both cases it works only for particular wavelength, so it should be specified while ordering. Laser radiation is delivered through FC/PCA fiber connector. It is focused onto PCA using lenses fixed inside the housing of THz emitter and THz detector.



Fig. 3. Femtosecond fiber laser LightWire FF50 with double output option



EKSPLA  
Savanoriu av. 231  
02300 Vilnius  
LITHUANIA

Ph.: +370 5 2649629  
Fax: +370 5 2641809  
sales@ekspla.com  
www.ekspla.com

Find local distributor at  
[www.ekspla.com](http://www.ekspla.com)

ISO 9001

CERTIFIED